

# Continuous VHU™ Perfusion for an order of magnitude increase in lentiviral vector production


M. Cattaneo<sup>1</sup>, A Rodenbrock<sup>2</sup>, RA Spanjaard<sup>1</sup> and S. Ansoorge<sup>2</sup> Affiliations: 1: Artemis Biosystems; 2: National Research Council Canada

## Introduction

- Lentiviral vectors (LV) represent a key tool for gene and cell therapy.
- The production of these vectors in sufficient quantities for clinical applications leads the field toward developing suspension processes more amenable to large-scale production.
- The HEK293 cell line employed grows in suspension, thus offering direct scalability, and produces a lentiviral vector in the 10<sup>6</sup> transduction units (TU)/mL range without optimization<sup>1</sup>.
- This study describes a LV production strategy using either a transient transfection or a stable inducible producer cell line grown in suspension using a novel perfusion process<sup>2</sup>.


## Methods

- A 3L bioreactor (BioBLU® 3c) was operated (37°C, pH 7.1, DO 40%, 80rpm) in perfusion mode using a novel cell retention device (VHU) (Figure 1).
- Cells were grown for 6 days with a perfusion rate of 0.5 volume vessel per day (vvd).
- Continuous VHU™ Perfusion increases the cell density up to 80x10<sup>6</sup> c/mL.
- Addition of inducers marked the beginning of the production phase and was considered Day 0. Continuous harvests were performed for up to 4 days post-induction (dpi) (Figure 2).



**Figure 1** The VHU™ Perfusion Bioreactor.


The integration of Viral Harvest Unit (VHU) (Top) and a low-shear pump (Quantum 600, Watson-Marlow) in TFF mode (Bottom) comprised the bioreactor perfusion system.



**Figure 2.** VHU™ Perfusion Process using a Stable Producer Cell Line


## Results

- VHU™ Process Intensification to increase HEK293 cell density up to 80 million cells/mL and with no cells in the perfusate (Figure 3).




**Figure 3.** VHU Perfusion compared to shake flask centrifugation. Viable cell density using VHU (red circles) versus shake flask mock perfusion (red triangles)

- Titer and Total Yield were 30-fold higher compared to Batch Mode (Figure 4).



**Figure 4.** (Left) Growth curve, Capacitance and Titer (GTA) for VHU™; (Right) Total LV titers for Batch vs. Perfusion with VHU™ Filter.

- VHU™ Perfusion provided a clarified stream of <16 NTU compared to <500 NTU for a batch culture (Figure 5 and Table 1).



**Figure 5.** Downstream Processing. LV Harvesting: Feed Flowrate: 3L/min; Permeate: 30mL/min; Clarified volume: 1.86 L; Functional titer: 7.58E+06 TU/mL

Process Mode	Turbidity (NTU)		
	Broth	Post-perfusion device	clarified
Batch	50-90	N/A	Target: ~5
Perfusion VHU	<500	<16	N/A

**Table 1.** Comparison of turbidity (NTU) between batch and perfusion cell retention devices such as the VHU™

## Conclusions

- The lentivirus was continuously harvested using the VHU™ for 4 days post induction.
- The VHU™ offers the potential to eliminate the clarification step.
- Batch vs. Perfusion: Perfusion ~3-times more expensive, 30-fold higher yield.
- Potential for additional higher yields using Process Intensification (80 million c/mL) in either transient transfection or stable producer cell line mode.

## REFERENCES

- Manceur A. et al, Scalable lentiviral vector production using stable HEK293SF producer cell lines. Hum Gene Ther Methods. 2017
- Cattaneo M. and Spanjaard R. Perfusion Filtration Systems, US Pat. 10,358,626 B2, 2019